

OPEN WIDE® Mouth Props

Safe oral access with superior comfort and stability

OPEN WIDE® Reusable Mouth Prop

Features

- ◆ Unique wrap around handle is designed to curve away from mouth.
- ◆ Doesn't block oral access and allows for full control of the prop and mandible at the same time.
- ◆ Cleverly designed "through the hole" receptor secures standard saliva ejector in place without the need for a second pair of hands.
- ◆ Made from highly crush-resistant hard plastic 'over moulded' with a soft, latex-free material that can be repeatedly autoclaved without getting gummy like other traditional wedges.



Available in pack of 3 colour coded props in small, medium and large sizes (paeds / adult)

OPEN WIDE® Disposable Mouth Prop

Single Patient Use
washable mouthprop perfect for home health, hospitals and dental clinics

Wedge-shape accommodates most patients

Rigid inner core adds stability and strength

Foam handle provides comfort and control

Handle-end can also be used for cheek retraction



Features

- ◆ Graduated wedge design with a padded handle for stability & comfort.
- ◆ Made from high density foam with a rigid core for added stability and strength making it difficult to bite through.
- ◆ For those with more aggressive bites a "thicker" version is available with 40% more foam for added strength.

Long favoured by Special Needs Dentists, the OWD is the preferred mouth prop to facilitate oral care in Nursing Homes, hospitals and home health settings.

The OWD is the perfect adjunct to the Surround® Toothbrush

